

HS Prime Taq DNA Polymerase

Product Name	Cat. No.	Size
HS Prime Taq DNA Polymerase	G-7000	250 Units X 1
	G-7001	250 Units X 2
	G-7002	250 Units X 4

Package Information

G-7000	1. HS Prime Taq DNA Polymerase (250 Units, 2.5 U/μl , 100μl X 1) 2. 10X Reaction buffer (with 20 mM MgCl ₂ , 1.0 ml X 1) 3. 10 mM dNTPs Mixture (2.5 mM each of dATP, dCTP, dGTP and dTTP. 0.5 ml X 1)
--------	--

Description

HS Prime Taq DNA Polymerase is designed for hot-start PCR, a technique that enhances the specificity, sensitivity and yield of DNA amplification. In addition, the enzyme provides the convenience of reaction set-up at room temperature. The enzyme is inactivated at room temperature, avoiding extension of non-specifically annealed primers or primer dimers and providing higher specificity of DNA amplification. The functional activity of the enzyme is restored during 10 minute incubation at 94°C. The activated enzyme maintains the same functionality as Taq DNA Polymerase.: it catalyzes 5'→3' synthesis of DNA, has no detectable 3'→5' proofreading exonuclease activity.

Buffers and Reagent

Storage Buffer

20 mM Tris-HCl (pH 8.0), 100 mM KCl, 0.5 mM EDTA, 0.1 mM DTT, 0.5% Tween 20, 0.5% Nonidet P-40 and 50% Glycerol

10X Reaction Buffer

Contains Tris-HCl (pH 9.0), PCR enhancers, (NH₄)₂SO₄ and 20 mM MgCl₂

10 mM dNTPs Mixture

2.5 mM each of dATP, dCTP, dGTP and dTTP

Applications

Hot-start PCR, RT-PCR, Amplification of low copy or high range size DNA targets, Real-time PCR, Multiplex PCR, T-vector cloning.

◆ Research Use Only

◆ Store at -20°C

Protocol

Optimal reaction conditions, such as reaction temperature & times, and amount of template DNA, may vary and must be individually determined.

The following 20μl reaction volume can be used for PCR.

1. Thaw the 10X Reaction buffer and 10 mM dNTPs Mixture.
2. Prepare a mastermix.

Component	Volume	Final conc.
Sterilized D.W.	add up to 20μl	-
10X Reaction Buffer	2.0μl	1X
10 mM dNTPs Mixture	0.5~2.0μl	0.25~1.0 mM
Upstream Primer (10 pmoles/μl)	0.2~2.0μl	0.1~1.0 pmoles
Downstream Primer (10 pmoles/μl)	0.2~2.0μl	0.1~1.0 pmoles
HS Prime Taq DNA Pol. (2.5 units/μl)	0.2~0.4μl	0.5~1.0 unit
Template DNA	Variable	10 fg~1μg

★Amount of template DNA

Bacteriophage λ, cosmid, plasmid DNA: 10 fg~300 ng

Genomic DNA: 100 ng~1μg

3. Mix the mastermix and dispense appropriate volumes into PCR tubes. Centrifuge the PCR tubes in a microcentrifuge for 10 seconds.
4. Perform PCR using your standard parameters (3-step cycling).

Step	Temp. & Time		Cycles
	Temp.	Time	
Initial denaturation	95°C	3~5 min.	1
Denaturation	95°C	30 sec.	25~35
Annealing	x°C	30 sec.	
Extension	72°C	30~60 sec.	
Final Extension	72°C	5 min.	1

★for PCR products longer than 3~4Kb, use an extension time of approximately 1 min, per Kb DNA.

5. Separate the PCR products by agarose gel electrophoresis and visualize with EtBr or any other means.

★A DNA fragment which is amplified by HS Prime Taq DNA Polymerase has A-overhang, and it enables you to do cloning by using T-vector.