

1 kb DNA Marker

Product Name	Cat. No.	Size
1 kb DNA Marker	M-3000	0.5 ml X 1
	M-3002	0.5 ml X 5

Specification


Concentration: 77 ng / μ l

Recommended loading: 2.0 ~ 5.0 μ l / 5 mm lane width

Number of bands: 10

Size range: 500 ~ 10,000 bp

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA,
 5% Glycerol, 0.005% Bromophenol Blue,
 0.005% Xylene Cyanol



4.0 μ l loading / 5 mm lane width
 1.0% TAE agarose gel with stained EtBr

● **Research Use Only**

● **Store at -20°C**